


GUNTHER WEISS SCIENTIFIC GLASSBLOWING CO., INC.

PROVIDING SOLUTIONS FOR LACRIMAL BYPASS SURGERY (CDCR) SINCE 1963

2024 Product Catalog



LOCATION:

5295 NE Elam Young Parkway
Suite 140
Hillsboro, OR 97124

EMAIL: info@guntherweiss.com

WEB: www.guntherweiss.com

PHONE: 503-644-4056

FAX: 503-644-0975

FDA: Registration #3005489143

FDA: Product Codes "HMX" and "HNL"


FDA: Regulation #886.4350

ISO 9001:2015 AGS-US111219-1

CONTENTS

Title Page.....	1
Contents.....	2
Part Number Guide.....	3
2024 Pricing Guide.....	4-5
10000 Jones Tube.....	6-7
20000 JT Gladstone-Putterman.....	8
50000 JT Callahan-Cox.....	9
60000 Specialty Jones Tubes.....	10-12
70000 Jones Tube Sets.....	13-14
80000 Weiss Glass Cannulas	15
90000 Accessories.....	16-17
How to Measure a Jones Tube.....	18
Terms and Conditions.....	19

PART NUMBER GUIDE



EXAMPLES:

- **10000 - 3512** Jones Tube Standard 3.5x12mm
- **10110 - 4020** Jones Tube Frosted with Angle 4.0x20mm
- **11000 - 4532** Jones Tube with Suture Hole 4.5x32mm

2024 PRICING GUIDE


Jones Tube Standard.....	10000-XXXX = \$55.00
Jones Tube no Bevel	10001-XXXX = \$53.00
Jones Tube with Angle	10010-XXXX = \$61.00
Jones Tube Frosted	10100-XXXX = \$90.00
Jones Tube Frosted with Angle	10110-XXXX = \$95.00
Jones Tube with Suture Hole	11000-XXXX = \$76.00
Jones Tube with Suture Hole, Angled	11010-XXXX = \$82.00
Jones Tube Frosted with Suture Hole	11100-XXXX = \$108.00
Jones Tube Frosted with Suture Hole, Angled.....	11110-XXXX = \$116.00
Jones Tube 130 Degree Angle	10002-XXXX = \$69.00
Jones Tube 130 Degree Angle with Suture Hole..	11002-XXXX = \$92.00
Jones Tube 130 Degree Angle Frosted.....	10102-XXXX = \$103.00
Jones Tube 130 Degree Angle Frosted with..... Suture Hole	11102-XXXX = \$111.00
JT Gladstone-Putterman	20000-XXXX = \$84.00
JT Gladstone-Putterman with Angle	20010-XXXX = \$89.00
JT Gladstone-Putterman Frosted	20100-XXXX = \$119.00
JT Gladstone-Putterman Frosted with Angle	20110-XXXX = \$123.00
JT Gladstone-Putterman with Suture Hole.....	21000-XXXX = \$105.00
JT Gladstone-Putterman with Suture Hole, Angled..	21010-XXXX = \$111.00
JT Gladstone-Putterman Frosted with Suture Hole..	21100-XXXX = \$139.00
JT Gladstone-Putterman Frosted with Suture..... Hole, Angled	21110-XXXX = \$144.00

JT Callahan-Cox	50000-XXXX = \$61.00
JT Callahan-Cox Frosted	50100-XXXX = \$92.00
JT Callahan-Cox with Suture Hole	51000-XXXX = \$81.00
JT Callahan-Cox Frosted with Suture Hole	51100-XXXX = \$116.00
Jones Tube Double Flange	60000-XXXX = \$108.00
Jones Tube Double Flange with Suture Hole.....	61000-XXXX = \$119.00
Jones Tube Double Flange Frosted.....	60100-XXXX = \$132.00
Jones Tube Double Flange Frosted with Suture... Hole	61100-XXXX = \$137.00
Jones Tube for CPAP.....	60001-XXXX = \$121.00
Jones Tube for Reduced Flow.....	60002-XXXX = \$115.00
Standard Set Jones Tubes	70002-0000 = \$792.00
Custom Set Jones Tubes.....	70002-0001 = \$792.00
Assorted Set Jones Tubes	70003-0000 = \$743.00
Special Set Jones Tubes	70004-0000 = \$838.00
For additional styles of sets, see pages 13-14	
Weiss Glass Cannula STRAIKO	80000-STRAIKO = \$115
Weiss Glass Cannula LWVI	80000-LWVI = \$184.00
Weiss Glass Cannula MIS.....	80000-MIS = \$189.00
Weiss Glass Cannula Besson.....	80000-Besson= \$120.00
Jones Gold Dilator #1 (1.7&1.9mm Arms).....	90000-0001 = \$70.00
Jones Gold Dilator #2 (2.2&2.4mm Arms).....	90000-0002 = \$70.00
Jones Gold Dilator #3 (2.7&2.9mm Arms).....	90000-0003 = \$77.00
Jones Tube Sterilization Tray (Blue Anodized)	90000-0005 = \$524.00
Jones Tube Measuring Scale.....	90000-0007 = \$64.00

10000- JONES TUBE (JT)

JT Standard

10000-XXXX
\$55.00




JT with Suture Hole

11000-XXXX
\$76.00



JT with Angle

10010-XXXX
\$61.00



JT Frosted

10100-XXXX
\$90.00



JT Frosted with Angle

10110-XXXX
\$95.00



JT no Bevel

10001-XXXX
\$53.00



EXAMPLES:

- 10000-3512 Jones Tube 3.5x12mm
- 10110-4020 Jones Tube Frosted with Angle 4.0x20mm
- 11000-4532 Jones Tube with Suture Hole 4.5x32mm

Flange Sizes:

3.0mm, 3.5mm, 4.0mm
4.5mm, 5.0mm, 5.5mm

Overall Length:

9.0mm-40.0mm
**additional sizes are available upon request*

10000- JONES TUBE (JT)

JT Frosted with Suture Hole

11100-XXXX
\$108.00



JT with Suture Hole, Angled

11010-XXXX
\$82.00



JT 130 Degree Angle with Suture Hole

11002-XXXX
\$92.00




JT Frosted with Suture Hole, Angled

11110-XXXX
\$116.00



JT 130 Degree Angle

10002-XXXX
\$69.00



JT 130 Degree Angle Frosted with Suture Hole

11102-XXXX
\$111.00



JT 130 Degree Angle Frosted

10102-XXXX
\$103.00



Flange Sizes:

3.0mm, 3.5mm, 4.0mm
4.5mm, 5.0mm, 5.5mm

Overall Length:

9.0mm-40.0mm


**additional sizes are available upon request*

EXAMPLES:

- 10002-4024 Jones Tube 130 Degree Angle 4.0x24mm
- 11002-3528 Jones Tube 130 Degree Angle with Suture Hole 3.5x28mm
- 11100-3012 Jones Tube Frosted with Suture Hole 3.0x12mm

20000- JT GLADSTONE- PUTTERMAN

**JT Gladstone-
Putterman**
20000-XXXX
\$84.00




**JT Gladstone-
Putterman
Frosted with
Angle**
20110-XXXX
\$123.00




**JT Gladstone-
Putterman
Frosted with
Suture Hole**
21100-XXXX
\$139.00



**JT Gladstone-
Putterman with
Angle**
20010-XXXX
\$89.00



**JT Gladstone-
Putterman with
Suture Hole**
21000-XXXX
\$105.00



**JT Gladstone-
Putterman
Frosted with
Suture Hole,
Angled**
21110-XXXX
\$144.00



**JT Gladstone-
Putterman
Frosted**
20100-XXXX
\$119.00



**JT Gladstone-
Putterman with
Suture Hole,
Angled**
21010-XXXX
\$111.00



Flange Sizes:

3.0mm, 3.5mm, 4.0mm
4.5mm, 5.0mm, 5.5mm

Overall Length:

9.0mm-40.0mm
**additional sizes are
available upon request*

EXAMPLES:

- 20000-3521 JT Gladstone-Putterman 3.5x21mm
- 21000-4012 JT Gladstone-Putterman with Suture Hole 4.0x12mm


***The Gladstone Putterman Jones Tubes are available in a larger overall size
(3mm tubing instead of 2.2mm tubing) Please call for part number and pricing**

50000- JT CALLAHAN- COX TUBES

JT Callahan-Cox

50000-XXXX


\$61.00



JT Callahan-Cox with Suture Hole

51000-XXXX

\$81.00



JT Callahan-Cox Frosted

50100-XXXX

\$92.00



JT Callahan-Cox Frosted with Suture Hole

51100-XXXX

\$116.00



Flange Sizes:

4.0 mm, 4.5mm, 5.0mm, 5.5mm

Overall Length:

9.0mm-40.0mm

**additional sizes are available upon request*

EXAMPLES:

- 50000-4518 JT Callahan-Cox 4.5x18mm
- 51000-4012 JT Callahan-Cox with Suture Hole 4.0x12mm
- 51100-5520 JT Callahan-Cox Frosted with Suture Hole 5.5x20mm


***The Callahan-Cox Jones Tubes are made with larger tubing
(3.0mm tubing instead of 2.2mm)**

60000- SPECIALTY JONES TUBES

Jones Tube Double Flange

60000-XXXX


\$108.00



Jones Tube Double Flange with Suture Hole

61000-XXXX

\$119.00



Jones Tube Double Flange Frosted

60100-XXXX


\$132.00



Jones Tube Double Flange Frosted with Suture Hole

61100-XXXX

\$137.00



Flange Sizes:

3.0mm, 3.5mm, 4.0mm, 4.5mm

Overall Length:

9.0mm-40.0mm

**additional sizes are available upon request*

EXAMPLES:


- 60000-4018 Jones Tube Double Flange 4.0x18mm
- 61000-3012 Jones Tube Double Flange Frosted 3.0x12mm

60000- SPECIALTY JONES TUBES CONT.

Jones Tube for CPAP

60001-XXXX

\$115.00



Jones Tube for CPAP with Suture Hole

61001-XXXX

\$126.00



Jones Tube for CPAP Frosted

60101-XXXX


\$131.00



Jones Tube for CPAP Angled

60011-XXXX

\$137.00



Flange Sizes:

3.0mm, 3.5mm, 4.0mm, 4.5mm

Overall Length:

9.0mm-40.0mm

**additional sizes are available upon
request*


EXAMPLES:

- All CPAP tubes can have combined attributes
- 60001-4017 Jones Tube for CPAP 4.0x17mm
- 61001-3511 Jones Tube for CPAP Frosted 3.5x11mm

60000- SPECIALTY JONES TUBES CONT.


Jones Tube for Reduced Flow

60002-XXXX
\$115.00



Jones Tube for Reduced Flow with Angle

60012-XXXX
\$122.00



Jones Tube for Reduced Flow with Suture Hole

61002-XXXX
\$121.00



Jones Tube for Reduced Flow Frosted

60102-XXXX
\$131.00



Flange Sizes:

3.0mm, 3.5mm, 4.0mm, 4.5mm

Overall Length:

9.0mm-40.0mm

**additional sizes are available upon request*

EXAMPLES:

- All reduced flow tubes can have combined attributes
- 61002-4520 Jones Tube for Reduced Flow with Suture Hole 4.5x20mm
- 61112-3512 Jones Tube for Reduced Flow Frosted, Angled and with Suture Hole 3.5x12

70000- JONES TUBE SETS

Assorted Set

70003-0000
\$743.00

Set includes the following parts:

- 1ea. 3.0x10mm
- 1ea. 3.0x12mm
- 1ea. 3.0x14mm
- 1ea. 3.0x16mm
- 1ea. 3.5x11mm
- 1ea. 3.5x13mm
- 1ea. 3.5x15mm
- 1ea. 3.5x17mm
- 1ea. 4.0x15mm
- 1ea. 3.5x14mm Angled
- 1ea. 3.5x15mm Angled
- 1ea. 4.0x14mm Angled
- 1ea. 4.0x15mm Angled
- 1ea. Gold Dilator #1
- 1ea. Gold Dilator #2

Standard Set

70002-0000
\$792.00

Set includes the following parts:

- 2ea. 3.0x14mm
- 2ea. 3.5x12mm
- 2ea. 3.5x13mm
- 2ea. 3.5x14mm
- 2ea. 3.5x15mm
- 2ea. 4.0x15mm
- 1ea. 4.0x16mm
- 1ea. 4.0x18mm
- 1ea. Gold Dilator #1
- 1ea. Gold Dilator #2

Tube Style	Assorted Set	Standard Set
Jones Tube Standard	70003-0000 \$743.00	70002-0000 \$792.00
Suture Hole	71003-0000 \$977.00	71002-0000 \$1,040.00
Frosted	70103-0000 \$1,170.00	70102-0000 \$1,252.00
Frosted with Suture Hole	71103-0000 \$1,426.00	71102-0000 \$1,527.00

70000- JONES TUBE SETS

Special Set

70004-0000
\$838.00

Set includes the following parts:

- 1 ea. 3.0x10mm
- 1 ea. 3.0x11mm
- 1 ea. 3.0x12mm
- 1 ea. 3.0x13mm
- 1 ea. 3.0x14mm
- 1 ea. 3.0x15mm
- 1 ea. 3.5x09mm
- 1 ea. 3.5x11mm
- 1 ea. 3.5x17mm
- 1 ea. 3.5x18mm
- 1 ea. 4.0x12mm
- 1 ea. 4.0x13mm
- 1 ea. 4.0x14mm
- 1 ea. 4.0x17mm
- 1 ea. 4.0x19mm
- 1 ea. Gold Dilators #1 and #2

Custom Set

70002-0001
\$792.00

Set includes the following parts:

- 14 tube sizes of your choice
- 1 ea. Gold Dilator #1
- 1 ea. Gold Dilator #2

Tube Style	Special Set	Custom Set
Jones Tube Standard	70004-0000 \$838.00	70002-0001 \$792.00
Suture Hole	71004-0000 \$1,103.00	71002-0001 \$1,040.00
Frosted	70104-0000 \$1,332.00	70102-0001 \$1,252.00
Frosted with Suture Hole	71104-0000 \$1,628.00	71102-0001 \$1,527.00

80000- WEISS GLASS CANNULAS

Weiss Glass Cannula- STRAIKO

80000-STRAIKO
\$115.00



Weiss Glass Cannula- MIS

80000-MIS
\$189.00



Weiss Glass Cannula- LWVI

80000-LWVI
\$184.00



Weiss Glass Cannula- Besson

80000-Besson
\$120.00



90000 ACCESSORIES

**Jones Gold
Dilator #1**
(1.7&1.9mm arms)
90000-0001
\$70.00



**Jones Gold
Dilator #2**
(2.2&2.4mm arms)
90000-0002
\$70.00




**Jones Gold
Dilator #3**
(2.7&2.9mm arms)
90000-0003
\$77.00



Jones Tube Measuring Scale

90000-0007
\$64.00



90000 ACCESSORIES

Jones Tube Sterilization Tray (Blue Anodized)

90000-0005

\$524.00



Measuring Scale
(reference only)

HOW TO MEASURE

A JONES TUBE

There are two ways to measure the overall length of a Jones Tube:


Measuring an Angled Jones Tube

- All Jones Tubes are measured with a caliper prior to shipping.
- When measuring a Jones Tube with angle the flange end is aligned flatly against one jaw of the caliper and the caliper is closed until it just touches the beveled end of the tube.
- The Jones Tube below is a 19mm Jones Tube with angle.



Using the Blue Anodized Sterilization Tray or Measuring Scale (90000-0005, 90000-0007)

- Always make sure that the flange end is touching the end of the measure scale marked with the zero.
- The beveled end of the tube should be rotated so that the very extreme beveled tip of the tube is aligned with the side of the scale that aligns best with the tip.
- The tolerance for the Jones Tube is +/- 0.3 MM.
- The Jones Tube in the photo below is 24mm



These three circular holes are used for reference to measure the top of the Jones Tube (flange). Calipers can also be used for sizes beyond this scale.

TERMS AND CONDITIONS

1. **Shipping:** FOB origin, typically UPS ground, unless otherwise specified within two weeks of order. No shipments via US mail due to losses and an inability to track shipments.
2. **Payment terms:** Net 30, 1.5% per month over 30 days is preferred. We do accept VISA or MasterCard for a 3% processing fee.
3. **Overseas Payment terms:** All overseas orders will need to be paid in advance with a credit card. No shipments will be sent until payment is received.
4. **Return policy:** Because we wish to ensure that our customers receive only newly manufactured Jones Tubes, we will not accept returns on these products. If you find you have ordered incorrect product please refuse our delivery prior to receiving your order from our shipper and before opening the box. Otherwise, we are unable to accommodate returns. Items other than Jones Tubes may be returned for credit if in new condition.
5. **Preferred carrier:** UPS
6. **Additional Charge:** Customers requesting special shipping instructions outside of Gunther Weiss Scientific's "preferred carrier" shall be assessed an additional handling charge for their request.

